

LiquiBand[®] Exceed[™]XS

We've drawn the line for an even, consistent wound coverage.
















NO CLOG TECHNOLOGY

- ▶ Tip will not clog, allowing consistent flow of adhesive
- ▶ 0.4mL of adhesive per applicator delivers approximately 7cm of coverage

PURE OCTYL FORMULATION



- ▶ Elliptical shape allows for narrow or wide application
- Porous felt tip allows drip-free priming and even, consistent application

EASY TO USE APPLICATOR

- ▶ Winged applicator allows for safe and easy activation
- ▶ Provides controlled delivery of adhesive





LiquiBand°Exceed[™]Range


Offers the improved attributes of precision and controlled delivery with added advantages of our innovative tip design.

LiquiBand®Exceed™XS → 7cm Yield → No clog tip → Ideal for small incisions, laparoscopic port sites

Key Benefits

- High Viscosity³
 - ▶ Provides optimal flexibility
- Controlled Expression of Adhesive⁴
 - ▶ Applicator offers control during expression of adhesive
- Microbial Barrier⁵
 - Acts as a barrier against gram-positive, gram-negative and fungal microbes. Proven effectiveness against S. aureus, P.aeruginosa, E.coli, Candida albicans and MRSA
- No Clog Formulation⁶
 - ▶ Allows users to work at their own pace
 - ► The adhesive will remain in-tact throughout the wound healing process
- Maintains wound closure throughout healing process⁹
 - ▶ The adhesive will stay on the skin for 5-10 days and sloughs off naturally





LiquiBand[®] Exceed[™]XS

Size	Qty Per Box	Product Code
0.4g	6	72014020

LiquiBand[®] Exceed[™]

Size	Qty Per Box	Product Code
0.8g	6	72014002



T 1-855-GLUES-77E info@liquiband.com

W liquiband.com